

**KERN**<sup>®</sup> **KERN & Sohn GmbH**

Ziegelei 1

72336 Balingen-Frommern

Germany

☎ +0049-[0]7433-9933-0

FAX +0049-[0]7433-9933-149

@ info@kern-sohn.com

# Operating manual KERN Interface Adapter with Cable - for Extension Box

## KERN YKUP-13

Typ TYKUP-13-A

Version 1.0

2021-09

GB



You will find the current version of these instructions also online under:

<https://www.kern-sohn.com/shop/de/DOWNLOADS/>

Under the column Operating instructions

TYKUP-13-A-BA-e-2110\_ExtensionBox



# KERN Interface Adapter with Cable

Version 1.0 2021-09

## Installation Instructions for Extension Box

### Contents

<b>1</b>	<b>General hints .....</b>	<b>2</b>
1.1	Installation.....	2

## 1 General hints

**Cable length:** 0,15 m

<b>i</b>	<ul style="list-style-type: none"><li>• Via the KERN Extension Box up to 3 interfaces of your choice can be connected.</li><li>• Cascading is not possible!</li><li>• <b>Only KERN KUP adapters may be connected to the 15-pol-sub-D-plug of the device!</b></li></ul>
----------	--

### 1.1 Installation

- Switch off the device
- Plug in the Extension Box at the 15-pol-sub-D-plug of the device
- Switch on device
  - The Extension Box is automatically recognised by the appliance. However, it may take up to one minute before the interfaces can be used.